

MP Biomedicals

FastPrep-24 5G™

Access DNA/RNA & Proteins from Any Sample Type



Lyse sample in 40 seconds or less!


High Yields

Flexible Sample Holders.

- Ambient to Cryogenic
- 2 mL to 50 mL capacity

FastPrep-24™ 5^G

No Sample is too tough for the FastPrep-24 5G!



MP Biomedicals - It's What's Inside That Counts!



NEW FEATURES	BENEFITS
Increased Power Settings to 10 M/S	Highest Power Available Release analytes of interest into protective buffer quicker. Lyse Tougher Samples
Time: 1-120 Seconds Cycles: 1-9 Cycle Pause: 1-300 seconds	Less User Intervention True Walk-Away Operation
Touchscreen User Interface	Intuitive and Easy to Operate
>70 Recommended Programs	Validated, for Quicker Optimizations,
Manual Program-N-Store (Save up to 12 custom programs)	Saves Time and Knowledge
Export Data Function	Archive Assays and Results
USB Port	Easy Software Upgrades
Sound Interface	Audible Program Run Status Alerts
Strobe Effect	Visual Optimization Tool

Lyse Efficiently:

- Tissue, Tumors, Bones
- Bacteria gram + or -
- Soil
- Yeast, Fungi, Spores
- Plants, Seeds, Roots
- Feces

Cat. No.	Description
116005500	FastPrep-24 5G Tissue Homogenizer

Optimize your research with
MP Lysing Matrix and FastDNA® Purification Kits
for the most intact DNA, RNA and Proteins
without compromising the integrity of your sample!

© 2015 Thermo Fisher Scientific Inc. All rights reserved.

Trademarks used are owned as indicated at www.fishersci.com/trademarks.

In the United States:

For customer service, call 1-800-766-7000

To fax an order, use 1-800-926-1166

To order online: www.fishersci.com

BN1230153

In Canada:

For customer service, call 1-800-234-7437

To fax an order, use 1-800-463-2996

To order online: www.fishersci.ca

